

City of Philadelphia Public Health Department
2230 Cottman Avenue Philadelphia, PA 19149
BinaxNOW™ COVID-19 Ag Card External Quality Control Log

Abbott recommends that external positive and negative controls be run:

- Once with each new shipment received
- Once for each untrained operator
- When required by local, state, and/or federal regulations, accrediting groups, or your lab's Quality Control procedures

If the expected external control results are not obtained, do not report patient results. Contact Technical Service at 800-257-9525.

KIT LOT NUMBER: _____

| DATE | BINAXNOW™ COVID-19 Ag CARD KIT LOT/EXP | POSITIVE CTRL LOT/EXP | NEGATIVE CTRL LOT/EXP | POSITIVE RESULT | NEGATIVE RESULT | CONTROL LINE PRESENT: Y/N | TESTER'S INITIALS | CORRECTIVE ACTION / COMMENTS |
|------|--|--------------------------|--------------------------|--------------------|--------------------|------------------------------|----------------------|---------------------------------|
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |

SUPERVISOR REVIEWED BY: _____ **DATE:** _____