

City of Philadelphia Public Health Department
2230 Cottman Avenue Philadelphia, PA 19149
BinaxNOW™ COVID-19 Ag Card Internal Controls and Patient Record

LOT NUMBER _____ **EXP. DATE** _____ **NAME OF FACILITY/SITE:** _____

Record the Date, Patient's Name, Patient Test Result, Internal Control Results and the tester's initials.

Positive Internal Control = The pink-to-purple line at the "Control" position can be considered an internal positive procedural control.

Negative Internal Control = The background color in the window should be light pink to white within 15 minutes.

DATE	PATIENT NAME	PATIENT ID NUMBER	PATIENT RESULTS	CONTROL LINE PRESENT?		INTERNAL CONTROL RESULTS		COMMENTS	TESTER'S INITIALS
				NO	YES	+	-		

SUPERVISOR REVIEWED BY: _____

DATE: _____